50012336

# Heraeus


INSTRUMENTS

D

GB

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Wärmeschränke T 6030, T 6060, T 6120, T 6200, T 6420, T 6760  
Umluft-Trockenschränke UT 6060, UT 6120, UT 6200, UT 6420, UT 6760  
Operating Instructions  
**Series 6000**  
Heating Ovens T 6030, T 6060, T 6120, T 6200, T 6420, T 6760  
Air-Circulation Drying Ovens UT 6060, UT 6120, UT 6200, UT 6420, UT 6760

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




## **Table of Contents**

<b>1. Introduction</b>	<b>3</b>
<b>2. Installation and Connections</b>	<b>4 - 7</b>
<b>3. Operating Elements</b>	<b>8 - 21</b>
Air-flap adjustment	
Temperature registration	
Temperature control	
Thermal protection	
Timers	
Lighting, speed regulator	
Master switch	
<b>4. Commissioning</b>	<b>22</b>
<b>5. Maintenance and Repairs</b>	<b>23 - 24</b>
<b>6. Appendix</b>	
Table 3	<b>25</b>
Classification of heating ovens acc. to DIN-Standard 12 880, Part 1	
Technical data	<b>26 - 28</b>
Dimensioned sketch	<b>29</b>
Controller adjustment	<b>30</b>
Eurotherm controller with interface RS 232	<b>31</b>
Eurotherm controller with interface RS 485	<b>32</b>
Pt 100 connection	<b>33</b>

## 1. Introduction

Your new **heating/air-circulation drying oven** is characterized by a high degree of economic efficiency and outstanding quality.

Please read the following instructions carefully so that you'll be in a position to operate your oven without any difficulty.

### Special Notes

**The heating/drying ovens meet the requirements of DIN-Standard 12 880 and VDE-Regulation 0700, Part 1.**

**Radio-interference suppression: acc. to VDE Regulation 0875, Part 1.**

**Type of protection: IP 20 acc. to DIN-Standard 40050.**

**The thermal safety class as per section "Thermal Protection" (Page 14) is of particular importance for the installation and operation of the oven.**

**In addition to the above, please observe all relevant local safety regulations.**

**The heating/drying ovens must not be used for drying processes or other heat treatments releasing combustible vapors which might form explosive mixtures with the air. Furthermore, the ovens are not designed for the heat treatment of dangerous dusts or fibrous materials.**

## **2. Installation and Connections**

**Ovens of thermal safety class 0 (see nameplate) must be installed so as to permit permanent attendance. They may not be stacked or built into existing laboratory furniture, cf. section "Thermal Protection" (Page 14).**

### **Transport**

Care should be taken not to lift the oven by the door or the door handle.

### **Connection to the Power Supply**

The ovens come with a power supply cable and a shock-proof plug (16 A).

For line fusing, a 16 A delay-action fuse is required.

Prior to connecting the unit to the line, make sure that the power supply ratings are those stated on the nameplate. The connection must be made in accordance with the relevant local safety regulations.

Rated voltage, rated frequency, rated power and current are listed under "Technical Data" and on the nameplate.

### **Place of installation**

*The oven must be placed on a level surface.*

### **Adjustable feet**

Its position is locked in if the adjustable feet are turned. This requires two wrenches (size 24).


**Articulated rollers**

Models T/UT 6420 and T/UT 6760 are fitted with 4 articulated rollers. The lever which releases the roller is situated above the locking lever. After the oven has been positioned, press the locking levers of the rollers.

**Stability  
(models 6420, 6760)**

In order to ensure optimum stability, turn the front rollers **forward** as soon as the oven is in position, and press the locking levers down (see fig. 1).

**Close the door(s) when the oven is being moved!**

**Working temperature**

The minimum working temperature for the T ovens is room temperature + 5°C, and room temperature + 20°C for the UT ovens.

**Avoid installation in drafty locations. The ambient temperature should not permanently exceed 35°C.**

The rear wall of the oven is fitted with 2 tubular sockets. When the air flap is set for fresh-air operation, cold air is drawn in through the lower socket and hot air blown out again through the upper one. With ovens of thermal safety class 0, and with units which do not feature air-flap adjustment, the upper socket (air-outlet socket) is sealed.

**Distance from walls**

**Distance between the unit and the walls (25 mm to the ceiling and side\*) should be sufficient so as to not obstruct the outgoing flow of air and to prevent damage from the hot air and vapors.**

**Exception\***

Models 6030, 6420 and 6760 require a distance of 200 mm from the ceiling.

Air-inlet pipe

An air inlet pipe can be connected to supply the oven with fresh air. The pipe should be dimensioned large enough to provide the desired exchange of air.

Air-outlet pipe


The vapors produced in the oven during the heating phase can be expelled through an air-outlet pipe.

**Vapors and gases detrimental to health and the environment must be carried off by all means.**

For this purpose, a temperature-resistant and corrosion-proof air-outlet pipe should be used which can be connected to the air-outlet socket.

When several ovens beside each other are connected to a central air-extraction pipe, we recommend the installation of a draft interrupter (see illustration).

Bench installation

Heating ovens without a fan (T 6060 and T 6120) can also be built into laboratory benches.

**UT ovens may not be built in.**

For safety reasons, the laboratory benches must be of non-combustible materials in accordance with DIN-Standard 4102. Installation is carried out as follows:

- \* Connect the heating oven to the line.
- \* Fit it underneath the bench, taking care not to damage the power-supply cable.

A wall spacer provides for clearance between the rear wall of the laboratory bench and the oven. The clearance all the way around the unit must be at least 25 mm.


## Stacking

### Accessories required I

Model:	6060	6120	6200
Stock no.:	50027660	50027661	50027661

T and UT models 6060, 6120 and 6200 can be stacked as follows:

A max. of 2 on top of each other and 2 stacks side-by-side;  
minimum distance from wall 150 mm;  
clearance between stacks: 25 mm.


The procedure is carried out in the following manner:

Remove the plastic rivets.


Put down the facing.

Secure the stacking elements and angular facing with head screws.

Place the second oven on top of the first one. The stacking elements fit into the hollow feet of the upper oven to secure it in place.

### 3. Operating Elements from page 8 - 21

3.1) Air-flap adjustment

3.2) Temperature registration

3.3) Temperature control

3.4) Upper/lower limit safeguard

3.5) Timers

3.6) Lighting, speed regulator

3.7) Master switch

### 3. Operating Elements

#### 3.1 Air-flap adjustment


open

closed

#### 3.2 Temperature registration


Temperature indicator  
Models 6030 - 6200 only

dial thermometer

40°C to 300°C

Temperature recorder

round-chart  
temperature recorder,  
measuring range:

0 to 300 °C

The round chart recorder works without ink or color ribbon. The curve is recorded as the sharp needle point presses periodically against the pressure sensitive record chart. But for the replacement of the record chart, this type of recording system is maintenance-free.

The record charts are graduated from 0 - 300 °C. Future orders require that the article number be stated.

#### For 50 Hz:

1-day record chart, art. no. : 50028271  
7-day record chart, art. no. : 50028272  
30-day record chart, art.no. : 50028273

#### For 60 Hz:

20h record chart, art. no. : 50029584  
140h record chart, art. no. : 50029585  
25-day record chart, art. no.: 50029586

Minimum quantity: 1 set (100 record charts).

#### Replacement of the record chart:

The plexiglass cover is detached by a slight counterclockwise rotation. After removal of the record chart, the following operating elements become visible:

#### Stroke adjustment:

Setscrew for adjustment of the stroke: clockwise rotation gives a heavier stroke.

### Feed selection

This slide is used to set the desired feed.

- 1 day / 7 days / 30 days (at a line frequency of 50 Hz).
- 20 h / 5.8 days / 25 days (at a line frequency of 60 Hz).

The pointer tip is visible in the narrow slot beneath the stroke-adjustment screw, and appears every 6 seconds to register the measuring point. Depending on the time of disconnection from the line, it is also possible for the pointer to stop in the recording position.


**In order to avoid damage to the pointer, if it is in the recording position, make sure that a record chart has been inserted before the plexiglass cover is closed.**

After the desired feed has been set, the new record chart is inserted and pushed onto the four guide pins in the center. Clockwise rotation of the driving shaft sets the time at the circumference of the record chart to the white mark at the edge of the instrument. A slight clockwise turn will lock on the plexiglass cover. (Note the inscription "ELMES 12" in the center.)

### 3.3 Temperature Control

Mechanical temp. controller  
"Heratherm"  
Models 6030 - 6200 only

Temperature range:  
40 - 250 °C

#### Setting:

Use either a coin or a screwdriver to release the locking device of the temperature controller.

**Note the marking !** 

Select the desired working temperature with the control knob.


Tighten the locking device again.

**Note the marking !** 

The set working temperature, once reached, is kept constant with the controller opening and closing the heating circuit as required. The white pilot lamp indicates that the heater is working.

If working space temperature and dial readout do not match, the setting must be corrected (see Appendix).

Electronic temperature controller "Digicon"

electronic controller w. digital temp. display for temp. range: 20 - 300 °C

**Setting:**

The electronic controller features digital display of the set and actual temperature.

Press the control button and select the desired working temperature. The set value appears in the digital display.

When the control button is released, the actual value appears.

The LED (top left) indicates that the heating system is on.

Digicon S with interface

The "Digicon S" temperature controller controls the temperature of the heating oven in accordance with the setpoint. The actual temperature is indicated digitally. When the setpoint button (blue) is pressed, the desired temperature replaces the actual value in the display.

The interface module is located above the "Digicon S" temperature controller (in models T/UT 6030, 6420 and 6760 it is next to it).

Operation with external program

Connect connecting cable (diode plug) to socket "W<sub>ext</sub> (0..20mA)".

0 mA = 0°C      20 mA = 300°C  
Admissible load R<sub>i</sub> = 500 Ω


This connection converts from the internal to the external setpoint program. There must be a bridge between pins 3 and 4 in the diode plug.

Measurement of the actual temperature

The actual temperature can be measured at the two "X<sub>ext</sub>." sockets as voltage of 0-10V.

0 Volt = 0°C      10 Volt = 300°C  
Admissible internal resistance R<sub>i</sub> = 1 kilo Ω

View of module with pin assignment



Connect only power sources / measuring amplifiers which are fed from a safety circuit (cf. VDE-Regulation 0100/410).

(Example: via safety transformer and VDE-Regulation 0551.)

Microprocessor-based temperature controllers

"Thermicon T"

"Thermicon P"  
(programmable)

For description and start-up, see enclosed operating instructions.

"Eurotherm  
type 818"

For description and start-up, see enclosed operating instructions.

### 3.4 Thermal protection

The requisite measures concerning thermal safety are stipulated in DIN-Standard 12 880, Part 1, Table 3 (see Appendix).

The thermal safety class which applies to these heating ovens is indicated on the nameplate.

#### No overtemperature protection

#### Thermal safety class 0:

No direct protection for the heating oven, its environment and load.

**Only admissible for non-dangerous loads.**

- If the temperature control system fails, the subsequent rise of temperature inside the oven will cause damage.
- Only admissible for supervised operation.
- The operation of these heating ovens must be checked at sufficiently short intervals.
- Heating ovens of this thermal safety class may not be built into laboratory furniture or stacked.

Supplementary equipment may be added to the heating ovens of safety class 0 if a temperature limit cut-out, temperature limit controller, or an upper/lower limit controller is installed.

## Temperature limit cut-out


### Thermal safety class 2

Protection of the heating oven, its environment and the load against inadmissible excess temperature.

**Only admissible for non-dangerous loads.**

### Function:

The temperature limit cut-out is operationally and electrically independent of the temperature control system, and disconnects the heater positively and from all poles.

### Adjustable temperature limit cut-out

When the control knob is set to the limit stop, the temperature limit cut-out functions as oven protection. If the temperature limit cut-out is set to a value somewhat higher than that selected on the controller, it functions as load protection.

Should the temperature limit cut-out trip the heating system (red warning lamp lights up):

1. Disconnect the oven from the line.
2. Investigate and eliminate the cause of the failure.
3. Press button to release temperature limit cut-out.
4. Turn the oven on again.

### Setting:

To ascertain the temperature at which the temperature limit cut-out operates, and thus to be able to adjust it properly, the oven must be put into operation. Set the control knob to the limit stop. After the selected working temperature has been reached, reset the temperature limit cut-out to the trigger position. If the cut-out comes into operation, the red warning lamp lights up, and the heating current is positively disconnected from all poles.

The optimum temperature limit cut-out adjustment is found if you turn back the control knob by at least two graduation marks. Then release the cut-out and finally set it approximately one graduation mark above the cut-off point. To release the lock, press the button thus unblocking the temperature limit cut-out. The red warning lamp goes out.

### Operational test:

The temperature limit cut-out should, at reasonable intervals, be tested for proper operation. It is recommended that this test be carried out routinely by the authorized operating personnel, e.g. prior to the start of a longer work process.

Connection for central monitoring system (potential-free contact)

Provision to connect customer-installed warning systems (e.g. telephone system, building control systems).

Function:


When the thermal safeguard (temperature limit cut-out, temperature limit controller) comes into action, or a timer switches the oven on or off, the appropriate message will be transmitted to the customer-installed warning system.

The potential-free contacts (1 make and 1 break contact) are designed for the following circuits:

- \* Power supply circuits:  
max. 250 V A.C., 3 A
- \* SELV/SELV-E circuits  
(cf. VDE-Regulation 0100/410):

25/50 V A.C., 4 A  
60/120 V D.C., 0.5 A

see also connection diagram

Temperature limit controller for excess temperature


Thermal safety class 3.1

Protection of the heating oven, its environment and the load against inadmissible excess temperature.

**Admissible only for non-dangerous loads.**

Function:

The temperature limit controller is operationally and electrically independent of the temperature control system. In case of failure, it takes over the control function.

Temperature limit controller

When the control knob is set to the limit stop, the temperature limit controller functions as oven protection. If the temperature limit controller is set to a value somewhat higher than that selected on the controller, it functions as load protection. Once the temperature limit controller begins to function and the red warning light comes on, proceed as follows:

1. Disconnect the oven from the line.
2. Investigate and eliminate the cause of the failure.
3. Turn the oven on again as described.

Setting:

To ascertain the temperature at which the temperature limit controller operates, and thus to be able to adjust it properly, the oven must be put into operation. Set the control knob to the limit stop (protection of the oven). After the selected operating temperature has been reached, reset the temperature limit controller to the switch point. If the controller comes into operation, the red warning lamp lights up.

The optimum setting of the temperature limit controller is found by turning back the control knob by approximately one graduation mark. The red warning lamp goes out.

Operational test:

The temperature limit controller should, at reasonable intervals, be tested for proper operation. It is recommended that this test be carried out routinely by the authorized operating personnel, e.g. prior to the start of a longer work process.

Upper/lower temperature limit controller  
functioning in case of excess and insufficient temperature 


### Thermal safety class 3.3

Protection of the heating oven, its environment and the load against inadmissible excess and insufficient temperature.

**Only admissible for non-dangerous loads.**

#### Function:

This controller complements the function of the temperature limit controller as per thermal safety class 3.1. With the flip switch, it can be turned on/off as desired and offers additional protection against insufficient load temperature if the controller fails.

Lower limit safeguard

Upper limit safeguard

Temperature limit controller with upper and lower limit safeguard.

#### Setting of the upper limit safeguard:

Set the desired temperature for take-over in the event of overheating as described under "temperature limit controller" (excess temperature).

#### Setting of the lower limit safeguard:

Turn on the white switch after the set temperature has been reached.

Pilot lamp in the switch lights up.

A slow turn of the knob will reveal the trigger position on the dial of the lower limit controller (the orange warning lamp lights up).

If you turn back the knob a little, the orange warning lamp will go out.

#### Operational test:

The upper/lower limit controller should, at reasonable intervals, be tested for proper operation. It is recommended that this test be carried out routinely by the authorized operating personnel, e.g. prior to the start of a longer work process.


**3.5 Timers:**  
Synchronous timer

**Function:**

Synchronous time switch  
for switching the oven on and off:

6 h duration at 50 Hz line frequency  
5 h duration at 60 Hz line frequency

24 h duration at 50 Hz line frequency  
20 h duration at 60 Hz line frequency

**24 h-Clock**

**Setting:**

Function of the mode switch

Position I: To switch on the oven after the set time has expired, or for operation without a synchronous timer.

Position II: To switch off the oven after the set time has expired.

**6 h-Clock**

For shut-off only.


Set the timer to the desired time by turning the knob to the right. Make sure to exceed the minimum adjustment angle of approx. 30°.

Shorter switching times are obtained if the timer is first turned to the right and then turned back to the desired time.

**Caution:**

The premounted dial (24 h) is valid for a line frequency of 50 Hz. For a line frequency of 60 Hz, the enclosed 20 h dial must be substituted.

The control button is released from the dial by impact from below (snug fit).


**Mounting of the 20 h dial**

- Set timer to 0 position
- Mount 20 h dial in 0 position
- Push on the control button.

**Daily program timer**

**Function:**  
To switch the oven on and off.

**Setting:**  
At the start of operation, set the actual time on the program dial by turning it to the right.


Select the desired switching times by locking down the riders.

The locked red dial segments indicate the switch-on times.

Minimum switching frequency 15 minutes.

**Weekly program timer**

**Function:**  
To switch the oven on and off.

Set the actual time and weekday on the program dial by turning it to the right.

1 = Mon., 2 = Tues., 3 = Wed., 4 = Thurs.,  
5 = Fri., 6 = Sat., 7 = Sun.

Select the desired switching time by locking down the riders.

Minimum switching frequency 2 hours.

**Mode converter:**

⌚ = automatic operation with timer.  
🔌 = manual operation: timer out of operation.

**Digital weekly program timer**


See separate operating instructions.

### 3.6 Switch combination:

Interior lighting  
with door window


**Caution!** Hot surface above door window can cause burns!

Interior lighting switch

Continuous operation not permissible!  
Maximum on time is 1 minute, followed by a minimum  
cool-down period of 9 minutes!

Speed regulation  
of air-circulation  
fan motor or fresh  
air blower

Infinitely variable from min. to max. speed. Function can  
also be switched off. For safety reasons, it cannot be  
switched off with types UT 6420/UT 6769 and fresh air  
blower-equipped units.

### 3.7 Master Switch:

Function:

Turn the knob to switch the oven on/off at all poles.

Turn on the master switch


- The green pilot lamp lights up.

## 4. Commissioning

### Shelves and supports

The shelves and their supports are secured inside the oven for transport. The transport guards must be removed before the oven is put into operation.

The supports can be inserted at any place in the perforated rows. The shelves slide over the supports in such a way that the fork-shaped stabilizer catches the lower side of the support.

**The working space proper is defined as that space which is 1/10 less than the overall chamber of the respective model. (cf. also DIN-Standard 12 880 Part 2).**

**Only this part of the chamber may be loaded.  
Do not load the bottom of the oven!**

### Note on loading

Don't arrange the load too closely on the shelves (load only 70% of the area), so that air circulation is not impaired and uniform heating is ensured.

Close the door.


Turn on the master switch.

- The green pilot lamp lights up.

Set the temperature controller to the desired working temperature.

After the desired temperature has been reached, set the temperature safeguard.

### Caution

The door surface is hot

- around the windows
- around the central edge of the doors (only with model 6760)

## 5. Maintenance and Repairs


The oven should, at reasonable intervals, be tested for proper operation. This applies especially to the temperature limit cut-out/controller.

**The safety of the oven and its conformity with the relevant standards will be guaranteed only if maintenance work, repairs and/or modifications are carried out by authorized personnel.**

Before any maintenance or repair work is carried out, the oven must be disconnected from the line.

### Readjusting the door

If the door no longer closes properly, adjust the bearing pin and the catches:

#### 1. Adjusting the catch

- Loosen the catch nut M6 with a wrench (size 10) and turn the catch clockwise for at least one complete rotation. After this, return the catch to the same position. Retighten nut M6.

#### 2. Adjusting the bearing pin

- Loosen the Phillips screw of the bearing pin. Use a small spindle to prevent the bearing pin from rotating also. After you've loosened the screw, turn the bearing pin with the spindle, and determine the desired setting for height and depth. Retighten the Phillips screw and check if the door fits tight.

### Replacing the gasket

Remove the defective gasket and install a replacement. After that it is imperative that the unit be tested for tightness. Make sure to reinstall the clamps.

### Removal of the internal fittings

Use a Phillips screwdriver to remove the internal fittings in the working chamber.

- Unscrew the bottom plate
- Unscrew the lateral cooling baffles
- \*Unscrew the rear wall covering
- Remove the individual internal fittings

\*only with UT ovens

Hints on care

Inner casing:

Commercially available detergents in small quantities should be used, but **no acids, no chloric solvents or saline solutions.**

Control elements (modules):

It is best to wipe these only with a damp cloth.

Spare parts

If you have complaints or wish to order for spare parts, please state the data on the nameplate.

A list of spare parts can be requested either from your local Heraeus representative or directly from the factory. State the model type and serial number.

Wiring diagrams

BR 6000 220 V/I A.C.	500 11 620
BR 6000 110 V/I A.C.	500 30 370
BR 6000 380 V/3 A.C.	500 11 621
BR 6000 220 V/3 A.C.	500 30 371

**6. Appendix - Table 3, from DIN-Standard 12 880**

Safety devices against malfunction of the temperature control circuit.

Classification of heating ovens according to DIN-Standard 12 880 Part 1, table 3.

The class describes the type and scope of protection, as well as the safety measures required to deal with a simple defect in the temperature control circuit.

Class	Purpose	Scope	Safety device acc. to Sub-section 5.5	Safety measures
0	-	-	-	Only attended operation*) with harmless materials is permitted; or overheating is excluded by constructional measures.
1	Protection of oven	In the event of a defect, no danger emanates from the oven.	Temperature limiter or thermal fuse	Special safety measures depending on the intended use.
2	Protection of oven, environment and material being treated.	In the event of a defect, no danger emanates from the oven or from the material being treated.	Adjustable temperature limit cut-out	
3		In the event of a defect, the material being treated is protected against overheating and/or undercooling (e.g. incubator)	Adjustable temperature limit controller, coming into action in case of: excess temperature insufficient temperature excess temperature and insufficient temperature	
3.1				
3.2				
3.3				

\*) If operation is supervised, the oven must be checked at reasonably short intervals to make sure it works correctly.

## Technical Data

	6030		6060		6120		6200		6420		6760	
	T	UT	T	UT	T	UT	T	UT	T	UT	T	UT
<b>Materials:</b>												
Inner casing	1.4301 *)		1.4301 *)		1.4301 *)		1.4301 *)		1.4301 *)		1.4301 *)	
Shelves	Chromium-plated steel		Chromium-plated steel		Chromium-plated steel		Chromium-plated steel		Chromium-plated steel		Chromium-plated steel	
Door gasket	Sil. rubber		Sil. rubber		Sil. rubber		Sil. rubber		Sil. rubber		Sil. rubber	
Heating	1.4435 *)		1.4435 *)		1.4435 *)		1.4435 *)		1.4435 *)		1.4435 *)	
Outer casing **)												
*) Stainless and acid-resistant steel												
**) Galvanized, sheet-steel outer coating with heat-resistant polyester varnish.												
<b>Colours:</b>												
Outer casing, grey-white	RAL 9002		RAL 9002		RAL 9002		RAL 9002		RAL 9002		RAL 9002	
Modules, agate-grey	RAL 7038		RAL 7038		RAL 7038		RAL 7038		RAL 7038		RAL 7038	
Operating elements, light-blue	RAL 5012		RAL 5012		RAL 5012		RAL 5012		RAL 5012		RAL 5012	
<b>Dimensions (mm)</b>												
see dimensioned sketch P.29												
<b>External dimensions</b>												
* Casing												
T <sub>A</sub> incl. door handle	610		610	647	610	647	790	827	790		790	
B <sub>A</sub>	552		744	744	895	895	895	895	744		1200	
H <sub>A</sub>	552		552	552	696	696	816	816	1813	1838	1813	1838
<b>Internal dimensions</b>												
T <sub>I</sub>	370		370	339	370	323	550	503	550	522	550	522
B <sub>I</sub>	352		403	403	554	554	554	554	544	544	1000	1000
H <sub>I</sub>	231		380	380	524	524	644	644	1366	1319	1366	1319
<b>Working-chamber volume (l)</b>												
	30		57	52	107	94	196	180	409	375	751	689

## Technical Data

	6030		6060		6120		6200		6420		6760	
	T		T	UT	T	UT	T	UT	T	UT	T	UT
<b>Weights (kg)</b>												
Empty weight	40		50	53	65	75	92	100	153	163	223	241
max. charge	50		50		50		75		75		150	
max. shelf load:												
surface	20		20		20		40		40		40	
spot	15		15		15		20		20		20	
<b>Electrical data:</b>												
Rated voltage (V)	220		220		220		220		380		380	
Rated freq. (Hz)	50 / 60		50 / 60		50 / 60		50 / 60		50 / 60		50 / 60	
Rated current (A)	3.4		6.5	6.8	9.3	10.3	11.1	12.2	6.5	6.8	8.2	8.6
Rated power (kW)	0.75		1.4	1.49	2.0	2.25	2.4	2.65	3.6	3.8	5.4	5.6
<b>Times: (minutes)</b>												
Heating-up time (air-flap closed) from 25°C to:												
300°C	75		35	55	60	60	65	75	70	80	80	85
150°C	30		10	20	20	25	30	23	35	30	55	30
70°C	15		10	10	10	15	20	7	40	10	55	12
Recovery time (door opened for 1 min) to 99% of the initial value at:												
300°C	7		6	6	9	11	10	10	18	6	8	11
150°C	6		4	3	4	5	5	2	11	6	5	6
70°C	2		4	1	4	1	5	2	6	1	3	1
<b>Air changes per h</b>												
Air flap open												
300°C	33		14	33	25	27	21	28	21	21	12	12
150°C	34		12	35	21	29	15	41	16	19	14	18
70°C	14		7	36	19	32	11	47	11.5	47	10	10





## Technical Data

	6030		6060		6120		6200		6420		6760	
	T		T	UT	T	UT	T	UT	T	UT	T	UT
<b>Max. fresh air quantity (m<sup>3</sup>/h)</b> (air flap open)												
300°C	1.4		1.06	2.5	3.41	3.68	5.04	6.7	10.3	10.3	10.3	10.3
150°C	1.45		0.91	2.65	2.86	3.96	3.6	9.8	7.8	9.8	12	15.5
70°C	0.6		0.53	2.73	2.6	4.36	2.64	11.27	5.6	9.8	8.6	8.6
with UT models at max. fan power												
quantity of recircu- lated air (at 25°C)												
max. fan power in m <sup>3</sup> /min	-		-	2.3	-	6.5	-	4.1	-	6.4	-	7
<b>*Holding power  (Wh/h)</b>												
air flap closed												
at 300°C	415		550	900	780	1280	1000	1570	1800	2500	2600	3300
150°C	150		200	410	265	580	350	695	610	1100	890	1250
70°C	55		75	185	85	265	120	325	210	420	280	520
<b>Temperatures**:</b>												
Rated temp.	300°C		300°C		300°C		300°C		300°C		300°C	
Variation in space (°C), air flap closed												
300°C	± 5.3		± 5	± 6	± 6	± 7	± 5.5	± 5.6	± 3.6	± 3.6	± 4.8	± 4.4
150°C	± 2.0		± 4	± 3	± 4	± 2.4	± 4	± 3.1	± 3.5	± 2.9	± 4.9	± 3
70°C	± 1.4		± 2	± 1	± 2	± 0.9	± 1.8	± 1.3	± 2.9	± 0.8	± 2.4	± 1
Fluctuation at rated temp. (°C)	≤ 0.5		≤ 0.5		≤ 0.5		≤ 0.5		≤ 0.5		≤ 0.5	
Accuracy of temp. display Telethermometer in % of full- scale reading	± 3		± 3		± 3		± 3		± 3		± 3	
Digital display in % of full- scale reading	± 1		± 1		± 1		± 1		± 1		± 1	

\* corresponds to the heat load of the installation area

\*\* measured in accordance with DIN-Standard 12 880, Part 2.

# Dimensioned Sketches

6030, 6420, 6760


6060, 6120, 6200

**Adjustment of the  
mechanical con-  
troller "Heratherm"**

If the temperature in the working space does not correspond to that on the dial of the mechanical controller, the dial can be readjusted.

To do this, the following steps are necessary:


1. Remove the blue button with the dial
2. Push the blue button out of the dial (cf. illustration below)
3. Reposition the dial on the module
4. Compare the temperature value with the dial value and readjust
5. Remount the blue button

**Eurotherm Controller with Interface RS 232**

**Assignment of the 25-pin D-Sub socket:**

Pin	1	Shield
	2	Rx ( reception )
	3	Tx ( transmission )
	4	( connected to 5 )
	5	( connected to 4 )
	6	( connected to 8 and 20 )
	7	COM ( common, ground)
	8	( connected to 6 and 20 )
	20	( connected to 6 and 8 )


----- All unlisted pins are not wired -----

The interface pin assignment is designed to permit direct connection to an IBM or IBM-compatible computer if the computer is equipped with a 25-pin plug. If the computer interface has a 9-pin plug, a commercially available adapter (9-pin socket and 25-pin plug) can be used.

Upon delivery of the unit, the controller address is set to 00, and the Baud rate to 9600. For other addresses or other Baud rates, the controller must be reconfigured (see operating instructions for the EUROTHERM controller).

## Eurotherm Controller with Interface RS 485

Assignment of the 25-pin D-Sub socket:

Pin	1	Shield
	2	( connected to 12 )
	3	Tx + ( transmission + )
	4	( connected to 5 )
	5	( connected to 4 )
	6	( connected to 8 and 20 )
	7	COM ( common, ground)
	8	( connected to 6 and 20 )
	12	Rx + ( reception + )
	13	Rx - ( reception - )
	16	Tx - ( transmission - )
	20	( connected to 6 and 8 )



----- All unlisted pins are not wired -----

Upon delivery of the unit, the controller address is set to 00, and the Baud rate to 9600. For other addresses or other Baud rates, the controller must be reconfigured (see operating instructions for the EURO THERM controller).



**Caution:** In the case of some RS 422/485 interface cards for personal computers, it may be necessary to change the **polarity** of the signals Rx and Tx. If the controller does not respond to the computer interface, connect Rx+ of the PC interface to Tx- of the controller, and Tx+ of the PC interface to Rx- of the controller.

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